Disp-X

Bottle-top Dispenser

RAININ
Pipetting 360°

Operating Manual

Please read carefully before use,
and follow all operating and safety
instructions!

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METTLER TOLEDO

Contents

Safety Instructions	
Application	
Materials	
Operating Exclusions	
Operating Limitations	
Components	
Preparing for use	
Connecting the discharge tube	
Connecting the filling tube	
Mounting instrument on bottle	
Priming the instrument	
Dispensing	
Setting the volume	
Connecting drying tube (optional)	
Cleaning	
Cleaning/replacing the filling valve	
Cleaning/replacing the discharge valve	
Autoclaving	
Preparation prior to the first autoclaving	

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Read carefully before use!

Safety Instructions

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual before operation.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves.
3. Observe all safety specifications provided by reagent manufacturers.
4. When dispensing inflammable media, make sure to avoid the buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

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Application

This instrument is designed for dispensing liquids, observing the following limits:

- 15 to 40 °C of instrument and reagent
- vapor pressure up to 500 mbar
- density up to 2.2 g/cm³
- kinematic viscosity up to 500 mm²/s
(dynamic viscosity [mPas] =
kinematic viscosity [mm²/s] x density [g/cm³])

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Materials

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials: Borosilicate glass, Al_2O_3 , ETFE, PFA, FEP, PTFE, platinum-iridium; PP (closure cap).

Operating Exclusions

Never use this instrument with

- liquids attacking ETFE, FEP, PFA, PTFE or Al_2O_3 (e.g., dissolved sodium azide)
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- liquids which are decomposed catalytically by platinum-iridium (e.g., H_2O_2)
- fuming acids
- trifluoroacetic acid
- tetrahydrofuran
- explosive liquids (e.g., carbon disulfide)
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument.

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Operating Limitations

Concentrated hydrochloric acid and concentrated nitric acid, chlorinated and fluorinated hydrocarbons, or chemical combinations which form deposits may make the piston difficult to move or may cause jamming (e.g., concentrated saline solutions, highly crystallizing solutions or highly concentrated alkaline solutions).

When dispensing inflammable media, make sure to avoid the buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

Warning!

If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Immediately stop dispensing and follow cleaning instructions (see p. 19) or contact the manufacturer.

Note:

Compatibility of the instrument for this special application (e.g., trace material analysis) must be checked by the user or contact the manufacturer.

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
Components

- A Piston seat
- B Cylinder sleeve
- C Volume adjustment
- D Valve locking ring
- E* Discharge valve
- F* Discharge tube
- G* Closure cap
- H Locking nut
- I* Filling tube
- J* Filling valve
- K Air vent cap
- L Valve block
- M* Glass cylinder with protective sleeve
- N* Piston
- O Mounting tool
- P Drying tube (optional)
- Q Seal

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Preparing instrument for use

Attention!

The adapters supplied with the instrument are made of polypropylene (PP), and can only be applied for media which do not attack PP.

Warning!

Wear protective clothing, eye protection and gloves. Follow all Safety Instructions and observe the Operating Exclusions and Limitations (see pages 8, 9).

Always wear protective gloves when touching the instrument and the bottle, especially when using dangerous liquids. Carry the mounted instrument only as shown in figure.

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
Connecting the discharge tube

Warning!

Only the slotted-type discharge tube supplied with the instrument allows correct operation. Do not use any other type of discharge tube. Never use damaged or deformed discharge tubes.


1. Push the discharge tube as far as possible into the discharge valve.
2. Firmly finger-tighten the locking nut. Check the discharge tube for a tight fit.

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


Connecting the filling tube

1. For instruments 0-25 ml and 0-50 ml, filling tube can be connected.
2. Cut off filling tube diagonal.
3. Push the filling tube into the valve block as far as possible – with the diagonal end down.

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Mounting instrument on bottle

The instrument can be screwed directly onto a 45 mm screwthread. For other bottles, use the accompanying adapters.

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Priming the instrument

Attention!

For small bottles use a bottle stand to prevent tipping over.

Warning!

Reagent may drip off the discharge tube and closure cap.

- Hold the discharge tube and carefully remove the closure cap.

- Slide the closure cap backwards on the tube support, away from the discharge tube opening.

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

1. To avoid splashes when priming hold discharge tube on inner wall of a suitable recipient.
2. Gently pull the piston up (1 to 2 cm), then press it down rapidly.
3. Repeat this procedure until most of the air bubbles in the glass cylinder and discharge tube have been removed. Pay particular attention to critical areas where bubbles may lodge (see black arrows).

Note:

Before using the instrument for the first time, ensure it is rinsed carefully or discard first few samples dispensed.

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Dispensing

Warning!

Follow all Safety Instructions (see page 5-7) and observe the Operating Exclusions and Limitations (see pages 8, 9). The discharge tube must always point away from the user or any other person. Never press down the piston when the closure cap is pushed on. Reagent may drip off the discharge tube and closure cap.

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

Volumen einstellen

- Hold an appropriate collecting vessel under the discharge tube opening.
- Pull the piston up slowly and evenly to the upper stop to draw reagent into the glass cylinder.
- Push the piston down gently and evenly to dispense the reagent into the vessel.

Attention!

After use, always leave the piston in the down position.

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
Connecting drying tube (optional)

Note:


The screwthreads of the bottle and/or the threaded adapter should be sealed with PTFE tape.

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Cleaning

The instrument must be cleaned in the following situations to assure correct operation:

- if the piston becomes difficult to move
- after dispensing liquids which form deposits
- before changing the reagent
- prior to long term storage
- should liquid have accumulated in the closure cap
- prior to sterilization
- before replacing the valves

Warning!

Follow all Safety Instructions (see pages 5-7). The glass cylinder, valves, filling and discharge tubes contain reagent. Point the valves and tube openings away from your body. Wear protective clothing, eye protection and appropriate hand protection.

Note:

Never change the pistons of the instruments!

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

1. Push on the closure cap.
2. Empty and unscrew the instrument from the bottle and lift the instrument so that the filling tube is out of the liquid but still in the bottle.
3. Carefully tap the filling tube against the bottle from the inside so that the reagent runs back into the bottle.

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4. Remove the closure cap, and dispense the remaining reagent into the bottle by repeated piston operation.
5. Mount the instrument onto another bottle filled with an appropriate cleaning solution.
6. Flush the instrument.
7. Empty the bottle, flush it and fill it with water. Repeat step 6 and 7.
8. Empty the instrument.

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
9. Unscrew the piston seat (A).
10. Carefully pull the piston (N) out of the glass cylinder (M) by pulling on the piston seat only.
11. Carefully remove deposits at the edge of the glass cylinder with a screw driver. Clean piston and cylinder with a bottle-brush.
12. Carefully reassemble the instrument in the reverse order. Flush instrument with distilled water.

Note:

Insert the piston straight into the cylinder, not at an angle.

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

Cleaning/replacing the filling valve

1. Follow all cleaning instructions (see pages 19-21).
2. Remove the filling tube.
3. Unscrew the filling valve using the mounting tool. Remove the valve with its sealing washer. Make sure the sealing washer does not remain in the valve block.
4. Rinse the filling valve in the cleaning solution, and clean with soft brush. If valve ball is stuck slightly press with a plastic pipette tip to loosen it.
5. First screw the cleaned/new filling valve together with its sealing washer into the valve block by hand and then tighten it firmly using the mounting tool.

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

Cleaning/replacing the discharge valve

1. Follow all cleaning instructions (see pages 19-21).
2. Unscrew the discharge tube.
3. Remove valve locking ring.
4. Unscrew the discharge valve using the mounting tool. Remove the valve with its sealing washer. Make sure the sealing washer does not remain in the valve block.
5. Rinse the discharge valve in the cleaning solution, and clean with soft brush. If valve ball is stuck slightly press with a plastic pipette tip to loosen it.
6. First screw the cleaned/new discharge valve together with its sealing washer into the valve block by hand and then tighten it firmly using the mounting tool.
7. Slide on the valve locking ring making sure the two prongs properly align with the notches.

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Autoclaving

This instrument is autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes, according to DIN EN 285. It is the user's responsibility to ensure effective autoclaving.

Preparation prior to the first autoclaving

- Unscrew filling valve (J) and discharge valve (E).
- Replace built-in PTFE-rings with the supplied elastic O-rings made of FKM.
- First screw in the filling valve (J) by hand and then discharge valve (E) and tighten them firmly using the mounting tool.

Note:

Observe the reduced chemical resistance of FKM!

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

1. Clean the instrument before sterilization (see pages 19-21).
2. Remove the discharge tube (F) with removed closure cap (G), filling tube (I), and drying tube (P), and put all parts on a towel. **Avoid contact with hot metal surfaces.**
3. Lay the instrument **with piston in down position** on the towel.
4. Autoclave all parts.

Note:

Allow to cool down in the closed autoclave to avoid deformation through too rapid cooling. After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. Do not re-assemble the instrument until it has cooled down to room temperature (cooling time approx. 2 hours).

Note:

Mounting tool can be autoclaved at 121 °C.

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Checking the Volume

The accuracy and coefficient of variation of the instrument are determined gravimetrically as follows:

- Set to the nominal volume.
- Dispense distilled H₂O.
- Weigh the dispensed quantity on an analytical balance.
- Calculate the dispensed volume taking the temperature into account.
- Perform at least 10 dispensing and weighing operations.
- Calculate the accuracy (A%) and the coefficient of variation (CV%) by means of the formula used in statistical quality control. The procedure is described e.g., in DIN EN ISO 8655-6. Observe the operating manual of the balance manufacturer and the corresponding standards.

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Calculations (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n} \quad \begin{array}{l} x_i = \text{results of weighings} \\ n = \text{number of weighings} \end{array}$$

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$$\mathbf{VK\%} = \frac{100 s}{\bar{V}}$$

Trouble Shooting

Problem	Possible cause	Corrective action
Piston difficult to move	<ul style="list-style-type: none"> - Formation of crystals, dirty 	<ul style="list-style-type: none"> - Stop dispensing immediately. Loosen piston with circular motion, but do not disassemble. Follow all cleaning instructions (➤ p. 19-21).
Filling not possible	<ul style="list-style-type: none"> - Volume adjusted to minimum setting - Sticking filling valve 	<ul style="list-style-type: none"> - Set to required volume (➤ page 16). - Clean the filling valve. If necessary replace the valve with sealing washer (➤ page 22).
Instrument does not dispense liquid	<ul style="list-style-type: none"> - Discharge tube closure cap pushed on - Discharge tube not firmly connected or damaged 	<ul style="list-style-type: none"> - Remove closure cap. - Push on the closure cap. Hold the piston in place, and push the original manufacturer's discharge tube on firmly (➤ page 12). Replace a deformed or damaged discharge tube (➤ Accessories and Spare Parts).
Air bubbles in the instrument	<ul style="list-style-type: none"> - Reagent with high vapor pressure has been drawn in too quickly - Seal not inserted, knurled locking nuts not firmly connected 	<ul style="list-style-type: none"> - Slowly draw in reagent. - Double check that seal (Q) is in place, and that the knurled locking nuts at the discharge tube are firmly seated and finger tight.

	<ul style="list-style-type: none"> - Air not expelled from the instrument - Filling tube not firmly connected or damaged - Valves dirty, not firmly installed or damaged 	<ul style="list-style-type: none"> - Priming the instrument (► page 14) - Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube. - Cleaning procedure (► page 19-21). Tighten the valves using the mounting tool. If necessary, replace the valves and sealing washers.
Dispensed volume is too low	<ul style="list-style-type: none"> - Discharge tube not firmly connected or damaged - Filling tube not firmly connected or damaged - Valves not firmly installed or damaged - Sticking filling valve 	<ul style="list-style-type: none"> - Push the original manufacturer's tube on firmly (► page 12). - Cleaning procedure (► page 19-21). Push the filling tube on firmly. If this does not suffice, cut off approx. 1 cm replace filling tube (► page 13). - Clean the instrument (► page 19-21). Tighten the valves using the mounting tool, if necessary replace the valves and sealing washers (► page 22) - Clean or replace discharge tube, if necessary.

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*Final test values (A = Accuracy, CV = Coefficient of variation) relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C) of instrument, ambience and dist. water, and with smooth operation.

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Items supplied:

operating manual, 1 discharge tube, 1 mounting tool,
1 filling tube and adapters of PP of the following sizes:

Nominal vol.	Adapter, PP	Filling tube
0.5 - 5	25 mm, 28 mm, 32 mm,	250 mm
1 - 10	38 mm, 40 mm	
2.5 - 25	32 mm, 38 mm, 40 mm	500 mm**
5 - 50		

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Accessories and Spare Parts

Discharge tube

FEP, complete with tube support, locking nut, and closure cap. Pack of 2.

for nominal volume/size/length		Cat. No./ MT Ordering #
5, 10 ml	bent, 80 mm	DX-DT10/ 17011747
25, 50 ml	bent, 120 mm	DX-DT50/ 17011748

Filling tubes

with adapter to connect 2 filling tubes, FEP. Pack of 2.

Length x outer diam., mm	Cat. No./ MT Ordering #
250 x 6.9	DX-FT250/ 17011749
335 x 6.9	DX-FT335/ 17011750
390 x 6.9	DX-FT390/ 17012004

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Nominal volume 5, 10 ml **Cat. No. / MT Ordering #**
DX-FV10 / 17011757

Nominal volume 25, 50 ml **Cat. No. / MT Ordering #**
DX-FV50 / 17011765

Discharge valve with sealing washer

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Return for Repair

Attention!

Transporting of hazardous materials without a permit is a violation of federal law.

RAININ will not accept instruments that are not appropriately cleaned and decontaminated.

Therefore, contact RAININ and obtain return authorization **before** sending your instrument for service.

Return the instrument, with the Return Authorization Number prominently displayed on the outside of the package to the address provided with the Return Authorization Number. Include an exact description of the type of malfunction and the media used.

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Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

Subject to technical modification without notice. Errors excepted.

Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

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Dispenser selection chart

Reagent	Disp-X
Acetaldehyde	+
Acetic acid (glacial), 100%	+
Acetic acid, 96%	+
Acetone	+
Acetonitrile	+
Acetylacetone	+
Acrylic acid	+
Acrylonitrile	+
Adipic acid	+
Allyl alcohol	+
Aluminium chloride	+
Amino acids	+
Ammonium chloride	+
Ammonium fluoride	+
Ammonia, 20%	+
Ammonium sulfate	+
n-Amyl acetate	+
Amyl alcohol (Pentanol)	+
Aniline	+
Barium chloride	+
Benzaldehyde	+
Chlorobenzene	+

Reagent	Disp-X
Chlorobutane	+
Chromic acid, 50%	+
Chromosulfuric acid	+
Copper sulfate	+
Cumene (Isopropyl benzene)	+
Cyclohexanone	+
Decane	+
1-Decanol	+
Dibenzylether	+
Dichlorobenzene	+
Diethanolamine	+
Diethylamine	+
1.2 Diethylbenzene	+
Diethylene glycol	+
Dimethyl sulfoxide (DMSO)	+
Dimethylaniline	+
Dimethylformamide (DMF)	+
Diphenyl ether	+
Ethanol	+
Ethanolamine	+
Ethyl acetate	+
Monochloroacetic acid	+

Reagent	Disp-X
Nitric acid, 30%	+
Nitrobenzene	+
Oleic acid	+
Oxalic acid	+
Perchloric acid	+
Petroleum	+
Phenol	+
Phenylethanol	+
Phenylhydrazine	+
Phosphoric acid, 85%	+
Phosphoric acid, 85% + Sulfuric acid, 98%, 1:1	+
Piperidine	+
Potassium chloride	+
Potassium dichromate	+
Potassium hydroxide	+
Potassium permanganate	+
Propionic acid	+
Propylene glycol (Propanediol)	+
Pyridine	+
Benzene (Benzol)	+

Benzoyl chloride	+
Benzyl alcohol	+
Benzylamine	+
Benzylchloride	+
Boric acid, 10%	+
Bromobenzene	+
Bromonaphthalene	+
Butanediol	+
1-Butanol	+
n-Butyl acetate	+
Butyl methyl ether	+
Butylamine	+
Butyric acid	+
Calcium carbonate	+
Calcium chloride	+
Calcium hydroxide	+
Calcium hypochlorite	+
Chloro naphthalene	+
Chloroacetaldehyde, 45%	+
Chloroacetic acid	+
Chloroacetone	+
Ethyl methyl ketone	+

Formaldehyde, 40%	+
Formamide	+
Glycerol	+
Glycol (Ethylene glycol)	+
Glycolic acid, 50%	+
Hexanoic acid	+
Hexanol	+
Hydriodic acid	+
Hydrochloric acid, 20%	+
Isoamyl alcohol	+
Isobutanol (Isobutyl alcohol)	+
Isopropanol (2-Propanol)	+
Isopropyl ether	+
Lactic acid	+
Methanol	+
Methoxybenzene	+
Methyl benzoate	+
Methyl butyl ether	+
Methyl formate	+
Methyl propyl ketone	+
Mineral oil (Engine oil)	+
Pyruvic acid	+

Salicyl acid	+
Salicylaldehyde	+
Scintillation fluid	+
Silver acetate	+
Silver nitrate	+
Sodium acetate	+
Sodium chloride	+
Sodium dichromate	+
Sodium fluoride	+
Sodium hydroxide, 30%	+
Sodium hypochlorite	+
Sulfuric acid, 98%	+
Tartaric acid	+
Tetramethylammonium hydroxide	+
Triethanolamine	+
Triethylene glycol	+
Urea	+
Zinc chloride, 10%	+
Zinc sulfate, 10%	+

Disp-X is not designed for use with Hydrofluoric acid (HF). Consult the "Operating Exclusions" and "Operating Limitations" sections of the Operating Manual (see pages 9, 10).

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact RAININ.

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